


Reagent Data Sheet

MV30F02

Versie nummer: 2

Versie datum: 29-05-2019

r-biopharm®**Reagent Data Sheet
Flumequine ELISA
ID: 5101FLUM**

Component ID	Description	SDS
5000PLSAR	Microtiter plate One sealed Microtiter plate (12 strips, 8 wells each)	-
5101FLUMCO	Flumequine conjugate Lyophilized	X
5101FLUMS16	Flumequine Standard 50 - 10 - 5 - 1 - 0,5 - 0,1 - 0 ng/ml	X
5101FLUMS	Flumequine concentrated Standard 1000 ng/ml	X
5101FLUMAB	Flumequine antibody Lyophilized	X
5000DBLZ	Dilution buffer Ready to use	X
5000RB	Rinsing buffer 20x concentrated	X
5000TMB	TMB substrate solution Ready to use	X
5000SS	Stop solution 0,5 M sulphuric acid	X

X SDS attached to this sheet