


Reagent Data Sheet

MV30F02

Versie nummer: 2

Versie datum: 29-05-2019

r-biopharm®**Reagent Data Sheet
Neomycin ELISA
ID: 5111NEO**

Component ID	Description	SDS
5000PLRAM	Microtiter plate One sealed Microtiter plate (12 strips, 8 wells each)	-
5111NEOCO	Neomycin conjugate Lyophilized	X
5111NEOS16	Neomycin Standard 20 - 10 - 5 - 2.5 - 1.25 - 0.625 - 0 ng/ml	X
5111NEOS	Neomycin concentrated Standard 1000 ng/ml	X
5111NEOAB	Anti-Neomycin antibody Lyophilized	X
5000MZZKF	Dilution buffer Ready to use	X
5000RB	Rinsing buffer 20x concentrated	X
5000TMB	TMB substrate solution Ready to use	X
5000SS	Stop solution 0,5 M sulphuric acid	X

X SDS attached to this sheet

Reagent Data Sheet

MV30F02

Versie nummer: 2

Versie datum: 29-05-2019

