


**Reagent Data Sheet**

MV30F02

Versie nummer: 2

Versie datum: 29-05-2019

r-biopharm®**Reagent Data Sheet  
Ivermectin ELISA  
5141IVER**

Component ID	Description	SDS
5000PLSAR	Microtiter plate One sealed Microtiter plate (12 strips, 8 wells each)	-
5141IVERCO	Ivermectin conjugate Lyophilized	X
5141IVERS16	Ivermectin Standard 50 - 25 - 10 - 5 - 2.5 - 1.25 - 0 ng/ml	X
5141IVERAB	Anti-Ivermectin antibody Lyophilized	X
50004xLZ	Dilution buffer 4x concentrated	X
5141IVEREB	Ivermectin Extraction Buffer Ready to use	X
5000RB	Rinsing buffer 20x concentrated	X
5000TMB	TMB substrate solution Ready to use	X
5000SS	Stop solution 0,5 M sulphuric acid	X

X SDS attached to this sheet