


Reagent Data Sheet

MV30F02

Versie nummer: 1

Versie datum: 17-10-2019

r-biopharm®**Reagent Data Sheet****Colistin ELISA****5151COL**

Component ID	Description	SDS
5151COLPL	Microtiter plate	-
	One sealed Microtiter plate (12 strips, 8 wells each)	
5151COLCO	Colistin conjugate	X
	100x concentrated	
5151COLS	Colistin Standard	X
	Lyophilized	
5091CAPDB	Dilution buffer	X
	4x concentrated	
5000RB	Rinsing buffer	X
	20x concentrated	
5000TMB	TMB substrate solution	X
	Ready to use	
5000SS	Stop solution	X
	0,5 M sulphuric acid	

X SDS attached to this sheet