


**Reagent Data Sheet**

MV30F02

Versie nummer: 1

Versie datum: 30-03-2020

r-biopharm®

<b>Reagent Data Sheet Plus Cow's Milk ELISA 5171MILK</b>
--

Component ID	Description	SDS
5171PLKCAS	Microtiter plate	-
	One sealed Microtiter plate (12 strips, 8 wells each)	
5171MILKCO	HRP labeled monoclonal antibody directed against GMP	X
	100x concentrated	
5171MILKS	K-Casein Standard	X
	Corresponding to 4% Cow's milk, Lyophilized	
5171DB	Dilution buffer	X
	Ready to use	
5000RB	Rinsing buffer	X
	20x concentrated	
5000TMB	TMB substrate solution	X
	Ready to use	
5000SS	Stop solution	X
	0,5 M sulphuric acid	

X SDS attached to this sheet