



## The fast, easy & robust solution – RIDA® QUICK mycotoxin rapid tests

Minimize the risks of undetected mycotoxin contaminations.

RIDA® QUICK mycotoxin tests – the perfect and accurate solution for on-site testing.

### Benefits



Short testing times



Easy to perform



Quantitative



The optimal tool for on-site testing  
of raw materials



Innovative test reading  
with smartphone

## RIDA® QUICK mycotoxin tests

Mycotoxin contaminations of food and feed are a growing problem all over the world. They account for billions of USD in economic losses worldwide. Make sure you get reliable results when testing for this naturally occurring hazardous toxins.

RIDA® QUICK quantitative mycotoxin tests are the optimal tool to get fast and reliable information about the

mycotoxin contamination of your cereal raw materials. They are fast and combine simple testing protocols with minimal equipment requirements. Mycotoxin testing with RIDA® QUICK products is easy and robust at the same time. A short training on the tests is sufficient to generate reliable quantitative results.

### General information

- **Technology:** Immunochromatographic test strips for quantitative results
- **Format:** 20 individually packed test strips

### Quantitative results

Quantitative test evaluation in combination with the innovative RIDA® SMART APP (Art. No. ZRSAM1000). RIDA® QUICK SCAN test reading is also supported.

#### Aqueous extraction

Product	Article number	Incubation time	Detection range RIDA® SMART APP
RIDA® QUICK Aflatoxin RQS ECO	R5209	5 min	2 - 300 µg/kg
RIDA® QUICK DON RQS ECO	R5911	3 min	250 - 50000 µg/kg
RIDA® QUICK Fumonisin RQS ECO	R5606	5 min	300 - 10000 µg/kg
RIDA® QUICK T-2 / HT-2 RQS ECO	R5304	5 min	50 - 10000 µg/kg



#### Ethanol/methanol extraction

Product	Article number	Incubation time	Detection range RIDA® SMART APP
RIDA® QUICK Aflatoxin RQS	R5208	3 min	2 - 300 µg/kg
RIDA® QUICK Zearalenon RQS	R5504	5 min	50 - 1000 µg/kg

For more information please contact your local R-Biopharm distributor or our International Sales office:

[sales@r-biopharm.de](mailto:sales@r-biopharm.de)

#### R-Biopharm contacts:

**International Sales**  
Phone: +49 (0) 61 51 - 81 02-0  
Fax: +49 (0) 61 51 - 81 02-40  
E-mail: [sales@r-biopharm.de](mailto:sales@r-biopharm.de)

**Order Department**  
Phone: +49 (0) 61 51 - 81 02-0  
Fax: +49 (0) 61 51 - 81 02-20  
E-mail: [orders@r-biopharm.de](mailto:orders@r-biopharm.de)