


Reagent Data Sheet

MV30F02

Versie nummer: 2

Versie datum: 29-05-2019

r-biopharm®**Reagent Data Sheet
Okadaic Acid ELISA
5191OKA**

Component ID	Description	SDS
5000PLSAR	Microtiter plate One sealed Microtiter plate (12 strips, 8 wells each)	-
5191OKACO	Okadaic acid conjugate 100x concentrated	X
5191OKAS16	Okadaic acid Standard 10 - 5 - 2 - 1 - 0.5 - 0.2 - 0 ng/ml	X
5191OKAAB	Okadaic acid antibody Lyophilized	X
500010xAN	Dilution buffer 10x concentrated	X
5000RB	Rinsing buffer 20x concentrated	X
5000TMB	TMB substrate solution Ready to use	X
5000SS	Stop solution 0,5 M sulphuric acid	X

X SDS attached to this sheet