


**Reagent Data Sheet**

MV30F02

Versie nummer: 1

Versie datum: 04-01-2021

r-biopharm®**Reagent Data Sheet****Florfenicol ELISA****ID: 5091FLORF**

<b>Component ID</b>	<b>Description</b>	<b>SDS</b>
5091FLORFPL	Microtiter plate	-
	One sealed Microtiter plate (12 strips, 8 wells each)	
5091FLORFCO	Florfenicol conjugate	X
	100x concentrated	
5091FLORFS16	Florfenicol Standard	X
	0,8 - 0,4 - 0,2 - 0,1 - 0,05 - 0,025 - 0 ng/ml	
5000CI	Extraction solution I	X
	Ready to use	
5000CII	Extraction solution II	X
	Ready to use	
500010xAN	Dilution buffer	X
	10x concentrated	
5000RB	Rinsing buffer	X
	20x concentrated	
5000TMB	TMB substrate solution	X
	Ready to use	
5000SS	Stop solution	X
	0.5 M sulphuric acid.	

X SDS attached to this sheet