


**Reagent Data Sheet**

MV30F02

Version nr. 2

Version date: 20-05-2019



<b>Reagent Data Sheet</b> <b>Methyltestosterone ELISA</b> <b>ID: 5081MTES</b>
---

Component ID	Description	SDS
5000PLSAR	Microtiter plate One sealed Microtiter plate (12 strips, 8 wells each)	-
5081MTESCO	Methyltestosterone conjugate 100x concentrated	X
5081MTESS16	Methyltestosterone Standard 1 - 0.5 - 0.25 - 0.125 - 0.063 - 0.031 - 0 ng/ml	X
5081MTESAB	Methyltestosterone antibody 100x concentrated	X
500010xAN	Dilution buffer 10x concentrated	X
5000RB	Rinsing buffer 20x concentrated	X
5000TMB	TMB substrate solution Ready to use	X
5000SS	Stop solution 0,5 M sulphuric acid	X

X SDS attached to this sheet