


Reagent Data Sheet

MV30F02

Versie nummer: 2

Versie datum: 29-05-2019

r-biopharm[®]

Reagent Data Sheet
19-Nortetosterone ELISA
ID: 5081NOR

Component ID	Description	SDS
5000PLSAR	Microtiter plate One sealed Microtiter plate (12 strips, 8 wells each)	-
5081NORCO	Nortestosterone conjugate Lyophilized	X
5081NORS16	Nortestosterone Standard 2 - 1 - 0.5 - 0.25 - 0.125 - 0.0625 - 0 ng/ml	X
5081NORAB	Anti-17 β -Nortestosterone antibody Lyophilized	X
5000DBMZ	Dilution buffer Ready to use	X
5000RB	Rinsing buffer 20x concentrated	X
5000TMB	TMB substrate solution Ready to use	X
5000SS	Stop solution 0,5 M sulphuric acid	X

X SDS attached to this sheet