

Safety data sheet

according to 1907/2006/EC, Article 31



Printing date 09.08.2021

Version number 1

Revision: 09.08.2021

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Trade name: Enzytec™ D-Isocitric Acid Reagent 1

Article number: E1222E11

1.2 Relevant identified uses of the substance or mixture and uses advised against

No further relevant information available.

Application of the substance / the preparation: In vitro

1.3 Details of the supplier of the safety data sheet

Manufacturer/Supplier:

R-Biopharm AG

An der neuen Bergstraße 17

D-64297 Darmstadt

Further information obtainable from: e-mail: info@r-biopharm.de

1.4 Emergency telephone number:

CHEMTREC (worldwide): +1 703-741-5970 / 1-800-424-9300

CHEMTREC (United Kingdom): 44-870-8200418 and 44-2038073798 (toll free)

For additional worldwide local emergency numbers refer to section 16.

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture

Classification according to Regulation (EC) No 1272/2008



health hazard

Repr. 1B H360D May damage the unborn child.



Skin Irrit. 2 H315 Causes skin irritation.

Eye Irrit. 2 H319 Causes serious eye irritation.

2.2 Label elements

Labelling according to Regulation (EC) No 1272/2008

The product is classified and labelled according to the CLP regulation.

Hazard pictograms GHS07, GHS08

Signal word Danger

Hazard-determining components of labelling:

imidazole

Hazard statements

H315 Causes skin irritation.

H319 Causes serious eye irritation.

H360D May damage the unborn child.

Precautionary statements

P280 Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P308+P313 IF exposed or concerned: Get medical advice/attention.

P332+P313 If skin irritation occurs: Get medical advice/attention.

P405 Store locked up.

(Contd. on page 2)

GB

Safety data sheet

according to 1907/2006/EC, Article 31



Printing date 09.08.2021

Version number 1

Revision: 09.08.2021

Trade name: Enzytec™ D-Isocitric Acid Reagent 1

P501

Dispose of contents/container in accordance with local/regional/national/international regulations.

(Contd. of page 1)

2.3 Other hazards**Results of PBT and vPvB assessment****PBT:** Not applicable.**vPvB:** Not applicable.

SECTION 3: Composition/information on ingredients

3.2 Chemical characterisation: Mixtures**Description:**

Mixture: consisting of the following components.

Mixture of substances with nonhazardous additions.

Dangerous components:

CAS: 288-32-4	imidazole	Repr. 1B, H360D	1-2.5%
EINECS: 206-019-2		Skin Corr. 1C, H314	
		Acute Tox. 4, H302	

Additional information: For the wording of the listed hazard phrases refer to section 16.

SECTION 4: First aid measures

4.1 Description of first aid measures**General information:** Immediately remove any clothing soiled by the product.**After inhalation:**

In case of unconsciousness place patient stably in side position for transportation.

After skin contact: Immediately wash with water and soap and rinse thoroughly.**After eye contact:**

Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.

After swallowing: If symptoms persist consult doctor.**4.2 Most important symptoms and effects, both acute and delayed**

No further relevant information available.

4.3 Indication of any immediate medical attention and special treatment needed

No further relevant information available.

SECTION 5: Firefighting measures

5.1 Extinguishing media**Suitable extinguishing agents:** Use fire extinguishing methods suitable to surrounding conditions.**5.2 Special hazards arising from the substance or mixture**

No further relevant information available.

5.3 Advice for firefighters**Protective equipment:** No special measures required.

SECTION 6: Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

Wear protective clothing.

6.2 Environmental precautions: Dilute with plenty of water.**6.3 Methods and material for containment and cleaning up:**

Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).

Dispose contaminated material as waste according to item 13.

(Contd. on page 3)

Safety data sheet

according to 1907/2006/EC, Article 31



Printing date 09.08.2021

Version number 1

Revision: 09.08.2021

Trade name: Enzytec™ D-Isocitric Acid Reagent 1

(Contd. of page 2)

6.4 Reference to other sections

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

SECTION 7: Handling and storage

7.1 Precautions for safe handling

Open and handle receptacle with care.

No special precautions are necessary if used correctly.

Information about fire - and explosion protection: Keep respiratory protective device available.

7.2 Conditions for safe storage, including any incompatibilities

Storage:

Requirements to be met by storerooms and receptacles: No special requirements.**Information about storage in one common storage facility:** Not required.**Further information about storage conditions:** Keep container tightly sealed.**7.3 Specific end use(s)** No further relevant information available.

SECTION 8: Exposure controls/personal protection

8.1 Control parameters

Additional information about design of technical facilities: No further data; see item 7.

Ingredients with limit values that require monitoring at the workplace:

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

Additional information: The lists valid during the making were used as basis.

8.2 Exposure controls

Personal protective equipment:

General protective and hygienic measures:

Keep away from foodstuffs, beverages and feed.

Immediately remove all soiled and contaminated clothing

Wash hands before breaks and at the end of work.

Store protective clothing separately.

Avoid contact with the eyes and skin.

Respiratory protection: Not required.

Protection of hands:



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

Material of gloves

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

(Contd. on page 4)

Safety data sheet

according to 1907/2006/EC, Article 31



Printing date 09.08.2021

Version number 1

Revision: 09.08.2021

Trade name: Enzytec™ D-Isocitric Acid Reagent 1

(Contd. of page 3)

For the permanent contact gloves made of the following materials are suitable:

Nitrile rubber, NBR

Recommended thickness of the material: ≥ 0.11 mmValue for the permeation: Level ≤ 480
As protection from splashes gloves made of the following materials are suitable:

Nitrile rubber, NBR

Recommended thickness of the material: ≥ 0.11 mmValue for the permeation: Level ≤ 480
Eye protection:


Tightly sealed goggles

SECTION 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties

General Information

Appearance:

Form:	Fluid
Colour:	Colourless
Odour:	Undistinguishable.
Odour threshold:	Not determined.

pH-value at 20 °C: 7.1

Change in condition

Melting point/freezing point:	0 °C
Initial boiling point and boiling range:	100 °C

Flash point: Not applicable.

Flammability (solid, gas): Not applicable.

Decomposition temperature: Not determined.

Auto-ignition temperature: Product is not selfigniting.

Explosive properties: Product does not present an explosion hazard.

Explosion limits:

Lower:	Not determined.
Upper:	Not determined.

Vapour pressure at 20 °C: 23 hPa

Density at 20 °C:	1 g/cm ³
Relative density	Not determined.
Vapour density	Not determined.
Evaporation rate	Not determined.

Solubility in / Miscibility with water:

Fully miscible.

Partition coefficient: n-octanol/water: Not determined.

Viscosity:

Dynamic at 20 °C:	0.952 mPas
Kinematic:	Not determined.

(Contd. on page 5)

GB

Safety data sheet

according to 1907/2006/EC, Article 31



Printing date 09.08.2021

Version number 1

Revision: 09.08.2021

Trade name: Enzytec™ D-Isocitric Acid Reagent 1

(Contd. of page 4)

Solvent content:**Water:** 97.7 %**Solids content:** 100.0 %**9.2 Other information**

No further relevant information available.

SECTION 10: Stability and reactivity**10.1 Reactivity** No further relevant information available.**10.2 Chemical stability****Thermal decomposition / conditions to be avoided:**

No decomposition if used according to specifications.

10.3 Possibility of hazardous reactions No dangerous reactions known.**10.4 Conditions to avoid** No further relevant information available.**10.5 Incompatible materials:** No further relevant information available.**10.6 Hazardous decomposition products:** No dangerous decomposition products known.**SECTION 11: Toxicological information****11.1 Information on toxicological effects****Acute toxicity** Based on available data, the classification criteria are not met.**LD/LC50 values relevant for classification:****288-32-4 imidazole**

Oral LD50 880 mg/kg (mouse)

Primary irritant effect:**Skin corrosion/irritation**

Causes skin irritation.

Serious eye damage/irritation

Causes serious eye irritation.

Respiratory or skin sensitisation Based on available data, the classification criteria are not met.**Additional toxicological information:****CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction)****Germ cell mutagenicity** Based on available data, the classification criteria are not met.**Carcinogenicity** Based on available data, the classification criteria are not met.**Reproductive toxicity**

May damage the unborn child.

STOT-single exposure Based on available data, the classification criteria are not met.**STOT-repeated exposure** Based on available data, the classification criteria are not met.**Aspiration hazard** Based on available data, the classification criteria are not met.**SECTION 12: Ecological information****12.1 Toxicity****Aquatic toxicity:** No further relevant information available.**12.2 Persistence and degradability** No further relevant information available.**12.3 Bioaccumulative potential** No further relevant information available.**12.4 Mobility in soil** No further relevant information available.**Additional ecological information:****General notes:** Not hazardous for water.**12.5 Results of PBT and vPvB assessment****PBT:** Not applicable.

(Contd. on page 6)

Safety data sheet

according to 1907/2006/EC, Article 31



Printing date 09.08.2021

Version number 1

Revision: 09.08.2021

Trade name: Enzytec™ D-Isocitric Acid Reagent 1

vPvB: Not applicable.

(Contd. of page 5)

12.6 Other adverse effects No further relevant information available.

SECTION 13: Disposal considerations

13.1 Waste treatment methods

Recommendation

Must not be disposed together with household garbage. Do not allow product to reach sewage system.

Disposal must be made according to official regulations.

Uncleaned packaging:
Recommendation: Disposal must be made according to official regulations.

Recommended cleansing agents: Water, if necessary together with cleansing agents.

SECTION 14: Transport information

14.1 UN-Number

ADR, IMDG, IATA

Void

14.2 UN proper shipping name

ADR, IMDG, IATA

Void

14.3 Transport hazard class(es)

 ADR, ADN, IMDG, IATA
Class

Void

14.4 Packing group

ADR, IMDG, IATA

Void

14.5 Environmental hazards:

Not applicable.

14.6 Special precautions for user

Not applicable.

14.7 Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable.

Transport/Additional information:

Not dangerous according to the above specifications.

UN "Model Regulation":

Void

SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture
Directive 2012/18/EU
Named dangerous substances - ANNEX I None of the ingredients is listed.

National regulations:
Waterhazard class: Generally not hazardous for water.

15.2 Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

GB

(Contd. on page 7)

Safety data sheet

according to 1907/2006/EC, Article 31



Printing date 09.08.2021

Version number 1

Revision: 09.08.2021

Trade name: Enzytec™ D-Isocitric Acid Reagent 1

(Contd. of page 6)

SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

Relevant phrases

H302 Harmful if swallowed.

H314 Causes severe skin burns and eye damage.

H360D May damage the unborn child.

Department issuing SDS: QS 0049-6151-8102-0

Contact: QS 0049-6151-8102-0

Worldwide local emergency numbers

Argentina Spanish Local (City) Buenos Aires Americas 54-1159839431
 Australia English Local (City) Sydney Oceania 61-290372994
 Austria German Local (City) Vienna Europe 43-13649237
 Belgium French Local (City) Brussels Europe 32-28083237
 Brazil Portuguese Local (City) Rio De Janeiro Americas 55-2139581449
 Brazil Portuguese Local (City) Sao Paulo Americas 55-1143491359
 Brazil - Toll Free Portuguese Toll free Mobile-Enabled Americas 0800 892 0479
 Bulgaria Bulgarian Local (City) Plovdiv Europe 359-32570104
 Cayman Islands English Local (National) Americas 345-749-8392
 Chile Spanish Local (City) Santiago Americas 56 2 2581 4934
 China Mandarin Local (National) Asia 4001-204937
 Colombia Spanish Toll Free Americas 01800-710-2151
 Costa Rica Spanish Local (National) Americas 506-40003869
 Croatia Croatian Local (City) Zagreb Europe 385-17776920
 Czech Republic Czech Local (City) Prague Europe 420-228880039
 Denmark Danish Local (National) Europe 45-69918573
 Dominican Republic Spanish Local (City) Santo Domingo Americas (829) 956-7588
 El Salvador Spanish Local (City) San Salvador Americas 503 2136 7633
 Estonia Estonian Local (National) Europe 372-6681294
 Finland Finnish Local (City) Helsinki Europe 358-942419014
 France French National Europe 33-975181407
 Germany German Local (City) Frankfurt Europe 49-69643508409 0800-181-7059
 Greece Greek Local (City) Athens Europe 30-2111768478
 Grenada English Local (City) St George Americas 1 (473) 230-0165
 Guinea French Americas 224 660 71 03 00
 Hong Kong Toll Free Asia 800-968-793
 Hungary Hungarian Local (City) Budapest Europe 36-18088425
 Iceland Icelandic Local (City) Reykjavik Europe 354 539 0655
 India Hindi, Bengali, English Asia 000-800-100-7141
 Indonesia Indonesian Toll Free Asia 001-803-017-9114
 Ireland English Local (City) Dublin Europe 353-19014670
 Israel Hebrew Local (City) Tel Aviv Asia 972-37630639
 Italy Italian Local (City) Milan Europe 39-0245557031 800-789-767
 Japan Japanese Local (City) Tokyo Asia 81-345209637
 Latvia Latvian Lettish Local (City) Riga Europe 371-66165504
 Lithuania Lithuanian Local (City) Vilnius Europe 370-52140238
 Luxembourg French and German Europe 352-20202416
 Macedonia Macedonian Local (City) Skopje Europe 389 2 551 7456
 Malaysia Malay Local (City) Kuala Lumpur Asia 60-39212579 41-800-815-308
 Mexico Spanish Americas 01-800-681-9531
 Netherlands Dutch Europe 31-858880596
 New Zealand English Local (City) Auckland Oceania 64-98010034
 Nigeria Hausa Local (City) Lagos Europe 234 1 227 8883
 Panama Spanish Local (National) Americas 507-8322475
 Peru Spanish Local (City) Lima Americas 51-17071295
 Philippines Tagalog Local (City) Manila Asia 1-800-1-116-1020

(Contd. on page 8)

Safety data sheet

according to 1907/2006/EC, Article 31



Printing date 09.08.2021

Version number 1

Revision: 09.08.2021

Trade name: Enzytec™ D-Isocitric Acid Reagent 1

(Contd. of page 7)

Poland Polish Local (City) Warsaw Europe 48-223988029
 Portugal Portuguese Local (National) Europe 351-308801773
 Romania Romanian Local (National) Europe 40-37-6300026
 Russia Russian Toll Free Mobile-Enabled Asia 8-800-100-6346
 Saudi Arabia Arabic and English CHEMTREC DID Only Asia 966-8111095861
 Singapore English and Mandarin Asia 65-31581349 or 800-101-2201
 Slovakia Slovak Local (City) Bratislava Europe 421-233057972
 Slovenia Slovenian Local (City) Ljubljana Europe 38-618888016
 South Africa English Toll Free Africa 0-800-983-611
 South Korea Korean Toll Free Asia 003-0813-2549 and 080-822-1374
 Spain Spanish Local (City) Barcelona Europe 34-931768545 and 900-868538
 Sweden Swedish Local (City) Stockholm Europe 46-852503403
 Switzerland German and French and Italian Zurich Europe 41-435082011
 Taiwan Mandarin Local (City) Taipei Asia 886-2-7741-4207 and 00801-14-8954
 Thailand Thai Toll Free Asia 01-800-13-203-9987
 Trinidad and Tobago English Local (National) Americas 1-868-224-5716
 Turkey Turkish Local (City) Istanbul Europe 90-212-7055340
 Ukraine Ukrainian Europe 380-947101374
 United Kingdom English Local (City) London Europe 44-870-8200418 and 44-2038073798

Abbreviations and acronyms:

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)
 IMDG: International Maritime Code for Dangerous Goods
 IATA: International Air Transport Association
 GHS: Globally Harmonised System of Classification and Labelling of Chemicals
 EINECS: European Inventory of Existing Commercial Chemical Substances
 ELINCS: European List of Notified Chemical Substances
 CAS: Chemical Abstracts Service (division of the American Chemical Society)
 LC50: Lethal concentration, 50 percent
 LD50: Lethal dose, 50 percent
 PBT: Persistent, Bioaccumulative and Toxic
 vPvB: very Persistent and very Bioaccumulative
 Acute Tox. 4: Acute toxicity – Category 4
 Skin Corr. 1C: Skin corrosion/irritation – Category 1C
 Skin Irrit. 2: Skin corrosion/irritation – Category 2
 Eye Irrit. 2: Serious eye damage/eye irritation – Category 2
 Repr. 1B: Reproductive toxicity – Category 1B

*** Data compared to the previous version altered.**

GB