


Reagent Data Sheet

MV30F02

Versie nummer: 2

Versie datum: 29-05-2019

r-biopharm®**Reagent Data Sheet
Ractopamine ELISA
ID: 5061RACT**

Component ID	Description	SDS
5061RACTPL	Microtiter plate One sealed Microtiter plate (12 strips, 8 wells each)	-
5061RACTCO	Ractopamine conjugate Lyophilized	X
5061RACTS16	Ractopamine Standard 2 - 1 - 0.5 - 0.25 - 0.125 - 0.0625 - 0 ng/ml	X
5061RACTS	Ractopamine concentrated Standard 100 ng/ml	X
5000MP	Milk powder	-
50004xK	Dilution buffer 4x concentrated	X
5000RB	Rinsing buffer 20x concentrated	X
5000TMB	TMB substrate solution Ready to use	X
5000SS	Stop solution 0,5 M sulphuric acid	X

X SDS attached to this sheet