


Reagent Data Sheet

MV30F02

Versie nummer: 2

Versie datum: 29-05-2019

r-biopharm®

**Reagent Data Sheet****β-Agonist ELISA****ID: 5061BAG**

Component ID	Description	SDS
5000PLSAR	Microtiter plate	-
	One sealed Microtiter plate (12 strips, 8 wells each)	
5061BAGCO	β-Agonist conjugate	X
	Lyophilized	
5061BAGS16	Salbutamol Standard	X
	2 - 1 - 0.5 - 0.25 - 0.125 - 0.0625 - 0 ng/ml	
5061BAGS	Salbutamol concentrated Standard	X
	100 ng/ml	
5061BAGAB	β-Agonist antibody	X
	Lyophilized	
5000DBHZ	Dilution buffer	X
	Ready to use	
5000RB	Rinsing buffer	X
	20x concentrated	
5000TMB	TMB substrate solution	X
	Ready to use	
5000SS	Stop solution	X
	0,5 M sulphuric acid	

X SDS attached to this sheet