


Reagent Data Sheet

MV30F02

Versie nummer: 2

Versie datum: 29-05-2019

r-biopharm®**Reagent Data Sheet**
β-Agonist Fast ELISA
ID: 5061BAGFc

Component ID	Description	SDS
5061BAGFcPL	Microtiter plate One sealed Microtiter plate (12 strips, 8 wells each)	-
5061BAGFcCO	BAGFc Conjugate 100x concentrated	X
5061BAGFcS16	BAGFc standaard 2 - 1 - 0,5 - 0,25 - 0,125 - 0,0625 - 0 ng/ml	X
5061BAGFcS	BAGFc concentrated Standard 100 ng/ml	X
50004xK	Dilution buffer 4x concentrated	X
5000RB	Rinsing buffer 20 times concentrated	X
5000TMB	TMB substrate solution Ready to use	X
5000SS	Stop solution 0.5 M sulphuric acid.	X

X SDS attached to this sheet