


Reagent Data Sheet

MV30F02

Versie nummer: 2

Versie datum: 29-05-2019

r-biopharm®**Reagent Data Sheet****AMAZ ELISA****ID: 5091AMAZ**

Component ID	Description	SDS
5091AMAZPL	Microtiter plate	-
	One sealed Microtiter plate (12 strips, 8 wells each)	
5091AMAZCO	AMAZ conjugate	X
	100x concentrated	
5091AMAZS16	AMAZ-NP Standard	X
	4 - 2 - 1 - 0,5 - 0,25 - 0,125 - 0 ng/ml	
5091AMAZS	AMAZ Standard	X
	AMAZ standard 100 ng/ml	
5000PBS15	Dilution buffer	X
	Ready to use	
5161MGSDDB	Sample dilution buffer	X
	Ready to use	
5000RB	Rinsing buffer	X
	20x concentrated	
5000TMB	TMB substrate solution	X
	Ready to use	
		X
5000SS	Stop solution	
	0.5 M sulphuric acid.	

X SDS attached to this sheet