


**Reagent Data Sheet**

MV30F02

Versie nummer: 2

Versie datum: 29-05-2019

r-biopharm®

**Reagent Data Sheet  
Chloramphenicol ELISA  
ID: 5091CAP**

Component ID	Description	SDS
5000PLSAR	Microtiter plate One sealed Microtiter plate (12 strips, 8 wells each)	-
5091CAPCO	Chloramphenicol conjugate Lyophilized	X
5091CAPS16	Chloramphenicol standaard 2 - 0,5 - 0,2 - 0,1 - 0,05 - 0,025 ng/ml	X
5091CAPS	Chloramphenicol standaard 100 ng/ml	X
5091CAPAB	Chloramphenicol antibody Lyophilized	X
5091CAPDB	Chloramphenicol sample dilution buffer 4x concentrated	X
5091CAPRCB	Chloramphenicol reconstitution buffer Ready to use	X
5000RB	Rinsing buffer 20x concentrated	X
5000TMB	TMB substrate solution Ready to use	X
5000SS	Stop solution 0,5 M sulphuric acid	X

X SDS attached to this sheet

## Reagent Data Sheet

MV30F02

Versie nummer: 2

Versie datum: 29-05-2019

r-biopharm®

