


**Reagent Data Sheet**

MV30F02

Versie nummer: 2

Versie datum: 29-05-2019



**Reagent Data Sheet**  
**Chloramphenicol Fast ELISA**  
**ID: 5091CAPF**

Component ID	Description	SDS
5000PLSAR	Microtiter plate One sealed Microtiter plate (12 strips, 8 wells each)	-
5091CAPFCO	Chloramphenicol Fast conjugate Lyophilized	X
5091CAPFS16	Chloramphenicol Fast standaard 4 - 1 - 0,4 - 0,2 - 0,1 - 0,04 ng/ml	X
5091CAPFS	Chloramphenicol Fast standaard 100 ng/ml	X
5091CAPFAB	Chloramphenicol Fast antibody Lyophilized	X
5091CAPDB	Chloramphenicol sample dilution buffer 4x concentrated	X
5091CAPRCB	Chloramphenicol reconstitution buffer Ready to use	X
5000RB	Rinsing buffer 20x concentrated	X
5000TMB	TMB substrate solution Ready to use	X
5000SS	Stop solution 0,5 M sulphuric acid	X

X SDS attached to this sheet