


**Reagent Data Sheet**

MV30F02

Versie nummer: 2

Versie datum: 29-05-2019

r-biopharm®

**Reagent Data Sheet  
Tetracycline ELISA  
ID: 5091TC**

Component ID	Description	SDS
5091TCPL	Microtiter plate	-
	One sealed Microtiter plate (12 strips, 8 wells each)	
5091TCCO	Tetracycline conjugate	X
	Lyophilized	
5091TCS	Tetracycline Standard	X
	2 ng/ml lyophilized	
5091TCS1000	Tetracycline concentrated standard	X
	1000 ng/ml lyophilized	
50004xK	Dilution buffer	X
	4x concentrated	
5091MTC	MTC Buffer	X
	Ready to use	
5000RB	Rinsing buffer	X
	20x concentrated	
5000TMB	TMB substrate solution	X
	Ready to use	
5000SS	Stop solution	X
	0,5 M sulphuric acid	

X SDS attached to this sheet