


Reagent Data Sheet

MV30F02

Versie nummer: 2

Versie datum: 29-05-2019

r-biopharm®

| |
|---|
| Reagent Data Sheet Sulphamethazine ELISA ID: 5101SUL |
|---|

| Component ID | Description | SDS |
|--------------|---|-----|
| 5000PLSAR | Microtiter plate | - |
| | One sealed Microtiter plate (12 strips, 8 wells each) | |
| 5101SULCO | Sulphamethazine conjugate | X |
| | Lyophilised | |
| 5101SULS | Sulphamethazine Standard | X |
| | 10 ng/ml lyophilised | |
| 5101SULAB | Anti-Sulphamethazine antibody | X |
| | Lyophilised | |
| 5000DBMZ | Dilution buffer | X |
| | Ready to use | |
| 5000RB | Rinsing buffer | X |
| | 20x concentrated | |
| 5000TMB | TMB substrate solution | X |
| | Ready to use | |
| 5000SS | Stop solution | X |
| | 0,5 M sulphuric acid | |

X SDS attached to this sheet