


**Reagent Data Sheet**

MV30F02

Versie nummer: 2

Versie datum: 29-05-2019

r-biopharm®

<b>Reagent Data Sheet</b> <b>Multi-Sulfonamides ELISA</b> <b>ID: 5101SULM</b>
---

Component ID	Description	SDS
5000PLRAM	Microtiter plate	-
	One sealed Microtiter plate (12 strips, 8 wells each)	
5101SULMCO	Multi-Sulfonamide conjugate	X
	100x concentrated	
5101SULMS16	Sulfamethazine Standard	X
	5 - 2,5 - 0,8 - 0,5 - 0,25 - 0,125 - 0 ng/ml	
5101SULMAB	Multi-Sulfonamide antibody	X
	Lyophilized	
5000PBS	Dilution buffer	X
	Ready to use	
5000RB	Rinsing buffer	X
	20x concentrated	
5000TMB	TMB substrate solution	X
	Ready to use	
5000SS	Stop solution	X
	0,5 M sulphuric acid	

X SDS attached to this sheet