


Reagent Data Sheet

MV30F02

Versie nummer: 2

Versie datum: 29-05-2019

r-biopharm®

Reagent Data Sheet Aflatoxin B1 Sensitive ELISA ID: 5121AFBS

Component ID	Description	SDS
5121AFBSPL	Microtiter plate	-
	One sealed Microtiter plate (12 strips, 8 wells each)	
5121AFBSCO	Aflatoxin B1 conjugate	X
	100x concentrated	
5121AFBSS	Aflatoxin B1 sensitive Standard	X
	100 pg/ml lyophilized	
5000PBS	Dilution buffer	X
	Ready to use	
5000RB	Rinsing buffer	X
	20x concentrated	
5000TMB	TMB substrate solution	X
	Ready to use	
5000SS	Stop solution	X
	0,5 M sulphuric acid	

X SDS attached to this sheet