


Reagent Data Sheet

MV30F02

Versie nummer: 2

Versie datum: 29-05-2019

r-biopharm[®]**Reagent Data Sheet
Aflatoxin Total ELISA
ID: 5121AFT**

Component ID	Description	SDS
5000PLRAM	Microtiter plate One sealed Microtiter plate (12 strips, 8 wells each)	-
5121AFTCO	Aflatoxin Total conjugate 100x concentrated	X
5121AFBS16	Aflatoxin B1 Standard 0.2 - 0.1 - 0.05 - 0.025 - 0.0125 - 0.0625 - 0 ng/ml	X
5121AFTAB	Anti-Aflatoxin Total antibody 100x concentrated	X
5000DBHZ	Dilution buffer Ready to use	X
5000RB	Rinsing buffer 20x concentrated	X
5000TMB	TMB substrate solution Ready to use	X
5000SS	Stop solution 0,5 M sulphuric acid	X

X SDS attached to this sheet