


Reagent Data Sheet

MV30F02

Versie nummer: 2

Versie datum: 29-05-2019

r-biopharm®

Reagent Data Sheet Deoxynivalenol ELISA 5121DON
--

Component ID	Description	SDS
5000PLSAR	Microtiter plate One sealed Microtiter plate (12 strips, 8 wells each)	-
5121DONCO	Deoxynivalenol conjugate Lyophilised	X
5121DONS16	Deoxynivalenol Standard 10 - 5 - 2.5 - 1.25 - 0.625 - 0.3125 - 0 ng/ml	X
5121DONAB	Anti-Deoxynivalenol antibody Lyophilised	X
5000MZZK	Dilution buffer Ready to use	X
5000RB	Rinsing buffer 20x concentrated	X
5000TMB	TMB substrate solution Ready to use	X
5000SS	Stop solution 0,5 M sulphuric acid	X

X SDS attached to this sheet