


**Reagent Data Sheet**

MV30F02

Versie nummer: 2

Versie datum: 29-05-2019

r-biopharm®**Reagent Data Sheet  
Bacitracin ELISA  
5151BAC**

Component ID	Description	SDS
5151BACPL	Microtiter plate	-
	One sealed Microtiter plate (12 strips, 8 wells each)	
5151BACCO	Bacitracin conjugate	X
	100x concentrated	
5151BACS16	Bacitracin Standard	X
	20 - 10 - 5 - 2.5 - 1.25 - 0.625 - 0 ng/ml	
5000PBS	Dilution buffer	X
	Ready to use	
5161MGSDB	Sample Dilution buffer	X
	Ready to use	
5000RB	Rinsing buffer	X
	20x concentrated	
5000TMB	TMB substrate solution	X
	Ready to use	
5000SS	Stop solution	X
	0,5 M sulphuric acid	

X SDS attached to this sheet