


**Reagent Data Sheet**

MV30F02

Versie nummer: 1

Versie datum: 23-09-2019

r-biopharm®**Reagent Data Sheet****Gluten-Tec®****ID: 5071GT**

Component ID	Description	SDS
5171GTPL	Microtiter plate	-
	One sealed Microtiter plate (12 strips, 8 wells each)	
5171GTCO	Gluten conjugate	X
	100x concentrated	
5171GTS16	Gluten Tec® Standaard	X
	5 - 2,5 - 1,25 - 0,625 - 0,313 - 0,156 - 0 ng/ml	
5091CAPDB	Dilution buffer	X
	4x concentrated	
5000RB	Rinsing buffer	X
	20x concentrated	
5000TMB	TMB substrate solution	X
	Ready to use	
5000SS	Stop solution	X
	0.5 M sulphuric acid.	

X SDS attached to this sheet