




# Purchasing options of the RIDA®SMART APP – including a list of supported smartphones


Art. No. ZRSAM1000

The RIDA®SMART APP is a **complete new evaluation technology for the quantification of RIDA®QUICK mycotoxins test strips (LFD)**. The software application RIDA®SMART APP is the only smartphone-based method for quantitative evaluation of LFD world-wide. The unique feature of the RIDA SMART APP is full connectivity to any cloud etc.

Customers can buy individual components or sets depending on their customer needs:

- For **onsite testing** e.g. during harvesting in the fields a validated android smartphone and the RIDA®SMART APP is used (option 1, 2 and 6).
- The RIDA®SMART BOX is used for **high throughput measurements in the lab**.





# Purchasing options of the RIDA®SMART APP – including a list of supported smartphones

Art. No. ZRSAM1000

R-Biopharm AG validates new smartphones regularly to include the most recent and widespread models. During validation of a smartphone the dynamic range of the camera at different light conditions is evaluated and adjusted. Smartphone stands to simplify your workflow are also available from R-Biopharm. As smartphones have different sizes, the stands are customized.

Updates of this flyer with the latest list of android smartphone can be found on our website at <https://food.r-biopharm.com/de/produkte/rida-smart-app/>.



The list of validated android smartphones is outdated very fast due to technology advances; therefore, smartphone models are often only available on the market for 12 months. Please note, that some smartphones are not available worldwide.

**Table 1:** Availability of **validated** android smartphones and stands from R-Biopharm AG

Validated smartphone models (android)	R-Biopharm smartphone availability Art. No.	R-Biopharm stand Art. No.
Huawei Nexus 6P	-	-
Google Pixel	-	-
Google Pixel XL	-	-
Google Pixel 2	-	-
Motorola Moto G6	-	-
Xiaomi Redmi Note 7	-	-
Xiaomi Redmi Note 9S	-	ZRSAS-REDMI9S
Xiaomi Redmi Note 10S	-	ZRSAS-REDMI10S
Xiaomi Redmi Note 12	ZRSAREDMI12	ZRSAS-REDMI12
Motorola G 52	-	ZRSAS-MotorolaG52
Samsung Galaxy A23	-	ZRSAS-GalaxyA23