

SureFood® ALLERGEN ID/QUANT **NEW** portfolio from January 2018

December 2017

SureFood® ALLERGEN ID and QUANT Kits are proven real-time PCR systems for the detection and quantification of food allergens according to EU regulation 1169/2011.

We are always trying to make things as easy as possible for our customer. That is why we decided to combine several real-time PCR systems (see table). The internal amplification control (IAC) was integrated in the reaction mix and the volume of the Taq polymerase was increased (0.7 µl instead of 0.1 µl per reaction), which simplifies the handling of the kits.

The new test kits of the ALLERGEN series contain following reagents: Reaction Mix (inclusive IAC), Taq Polymerase, Positive Control, Standard DNA and Dilution Buffer. In this way the test systems can be used all in one for qualitative and quantitative analysis.

Starting from **1st January 2018** the following changes will be made:

Until December 2017 (+ 6 transitional period)		Starting from January 2018
SureFood® ALLERGEN Soya (S3101) SureFood® ALLERGEN QUANT Soya (S3201)	⇒	SureFood® ALLERGEN Soya (S3601)
SureFood® ALLERGEN Hazelnut (S3102) SureFood® ALLERGEN QUANT Hazelnut (S3202)	⇒	SureFood® ALLERGEN Hazelnut (S3602)
SureFood® ALLERGEN ID Peanut (S3103) SureFood® ALLERGEN QUANT Peanut (S3203)	⇒	SureFood® ALLERGEN Peanut (S3603)
SureFood® ALLERGEN ID Almond (S3104)	⇒	SureFood® ALLERGEN Almond (S3604)
SureFood® ALLERGEN ID Celery (S3105) SureFood® ALLERGEN QUANT Celery (S3205)	⇒	SureFood® ALLERGEN Celery (S3605)
SureFood® ALLERGEN ID Gluten (S3106) SureFood® ALLERGEN QUANT Gluten (S3206)	⇒	SureFood® ALLERGEN Gluten (S3606)
SureFood® ALLERGEN ID Walnut (S3107) SureFood® ALLERGEN QUANT Walnut (S3207)	⇒	SureFood® ALLERGEN Walnut (S3607)
SureFood® ALLERGEN ID Sesame (S3108) SureFood® ALLERGEN QUANT Sesame (S3208)	⇒	SureFood® ALLERGEN Sesame (S3608)
SureFood® ALLERGEN ID Mustard (S3109)	⇒	SureFood® ALLERGEN Mustard (S3609)
SureFood® ALLERGEN ID Fish (S3110)	⇒	SureFood® ALLERGEN Fish (S3610)
SureFood® ALLERGEN ID Lupin (S3111) SureFood® ALLERGEN QUANT Lupin (S3211)	⇒	SureFood® ALLERGEN Lupin (S3611)
SureFood® ALLERGEN ID Crustaceans (S3112)	⇒	SureFood® ALLERGEN Crustaceans (S3612)
SureFood® ALLERGEN ID Molluscs (S3113)	⇒	SureFood® ALLERGEN Molluscs (S3613)
SureFood® ALLERGEN ID Pistachio (S3114) SureFood® ALLERGEN QUANT Pistachio (S3214)	⇒	SureFood® ALLERGEN Pistachio (S3614)
SureFood® ALLERGEN ID Cashew (S3115)	⇒	SureFood® ALLERGEN Cashew (S3615)
SureFood® ALLERGEN ID Macadamia (S3116)	⇒	SureFood® ALLERGEN Macadamia (S3616)
SureFood® ALLERGEN ID Brazil Nut (S3117)	⇒	SureFood® ALLERGEN Brazil Nut (S3617)
SureFood® ALLERGEN ID Pecan (S3118)	⇒	SureFood® ALLERGEN Pecan (S3618)

remain available

- SureFood® ALLERGEN 4plex Soya/Celery/Mustard + IAC (S3401)
- SureFood® ALLERGEN 4plex Peanut/Hazelnut/Walnut+IAC (S3402)
- SureFood® QUANTARD Allergen 40 (S3301)

Thank you for your understanding. Please do not hesitate to contact us for any further questions.

With best regards,

CONGEN Biotechnologie GmbH
i.A. E. Lorenzen